



## Opmaster 506P

### A purpose-designed ophthalmic motorised table:

- Offers the widest range of height adjustment in any mobile operating table available today.
- Special ophthalmic features and 4-section flexibility make an ideal solution for busy day surgery units.
- Head section with articulating/levelling.
- Versatile dedicated mobile ophthalmic operating table.
- Height adjustment from 578mm to 978mm.

### Especially for ophthalmic procedures

- Motorised height range allows adjustment to create the ideal working position.
- Extended back section and choice of head sections give easy access to the patient while allowing maximum leg-room and space for microscope pedals.
- Motorised fine adjustment head section.
- Motorised trendelenburg.
- Motorised reverse trendelenburg and back section.
- Making patient positioning easy and strain-free for medical staff.

### Simplifies nursing duties

- Optional removable cot side rails give easy all round access to the patient.
- Convenient I/V pole holders remove the need for clamps.

### Cuts the strain of patient handling

- Motorised 4-section flexibility.
- Can be used as an operating table, transfer trolley, recovery bed and chair.
- Low elevation is a major benefit removing the strain of transfer and allowing the patient to sit directly onto the table when it is in the chair position.
- Easy roll castors provide maximum manoeuvrability.
- Balanced design and sturdy construction offer exceptional stability in every configuration.

### Enhanced patient comfort

- Has 75mm deep injection moulded, anti-static, easy-to-clean mattress for patient comfort.
- Mattress width allows patients to lie with their arms by their sides without the need for additional arm boards.
- Sculptured Rubens type head cushion effectively and comfortably restrains the head during procedures.

### Exceptional versatility

- Most versatile dedicated mobile ophthalmic operating table available today.
- Accessory bars and a full range of accessories including the Chang wrist support, sculptured head cushion and oxygen delivery system.

### Power when you need it

- Enclosed fully charged battery, in good condition will give a minimum of 1 hour continuous use.
- Simply recharge the battery in situ by direct mains, or use optional wall-mounted charging station.

### Dimensions

- Max height to mattress top: 978mm
- Min height to mattress top: 578mm
- Length: 1990mm - Length in chair position: 1250mm
- Width over mattress: 600mm
- Width over side rails: 750mm
- Mattress: 75mm - Total weight: 160kg

### Adjustments

- Trendelenburg: 20 degrees
- Reverse trendelenburg: 15 degrees
- Foot section: 87 degrees
- Head section: +36/-36 degrees
- Chair position: 60 degrees - Weight capacity: 225kg